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Stand der Information: 05/2011

Druckdatum: 11.11.2011

Polya BioStar PP Instructions for use

CE 0124 Dental disc on a thermoplastic acrylic basis

Polya BioStar PP is a thermoplastic acrylic milling blank designed for making long-lasting temporary crown and bridge work. Polya BioStar PP can be used in all CAD-CAM machines.

- very high breaking resistance
- no toxic agents
- free of residual monomer (below 0.3 %)
- adding of colour shades are no problem
- plaque resistant – high surface density
- benzyol peroxide free
- high translucence property
- perfect for telescopic crown

Composition:	Methyl methacrylate-based acrylic polymers	
Tensile strength	80 MPa	According to DIN EN ISO 527
Elastic limit	60 MPa	According to DIN EN ISO 527
Elongation at break	6 %	According to DIN EN ISO 527
Modulus of elasticity	3400 MPa	According to DIN EN ISO 527
Bending strength	106,0 MPa	According to DIN EN ISO 178
Water absorption	<32 µg/mm ²	Comparable with EN ISO 10477:1996

General: To process Polya BioStar PP, use only 2.5 mm tungsten carbide burs and 1.00 mm tungsten carbide burs. Please note that the improper use of milling tools, polishing brushes, steam jets and water baths can cause the material to overheat. This can result in damage to the material.

Notes on design: For crown and bridge work the following minimum values must be adhered to:

- Minimum wall thickness (cervical): 0,6 mm
- Minimum wall thickness (occlusal): 1,2 mm
- Cross section of connectors in the anterior region: 10 mm²
- Cross section of connectors in the posterior region: 16 mm²

In the case of posteriors, no more than two pontics should be placed between any two abutment teeth.

Notes on milling and finishing: Use a fine-toothed tungsten carbide bur to separate the milled framework from the blank. Cut the retaining bars through carefully and without applying pressure. The Polya BioStar PP restorations can be prepolished with suitable silicone polishes and goat hair brushes. To achieve a high shine polish, use the diamond polishing pastes (Polirapid P3) recommended for the purpose. Thorough polishing is absolutely necessary to avoid plaque build-up.

Characterisation: Composites can be used to individually characterize the aesthetics of a restoration. To enable this characterisation, the corresponding areas are reduced to allow a smooth transition to the base material. In doing so, not more than 0,3 mm – 0,5 mm of material should be removed. Only fine-toothed tungsten carbide burs should be used for grinding. The ground surface must be carefully cleaned and wetted with a suitable bonder. For details of handling the veneering materials consult the manufacturer's instructions.

Cleaning: With all plastics and composites used for making frameworks or applying veneers, steam cleaning can cause undue heat and pressure and should be avoided as a matter of course. Alternative methods of cleaning are physical cleaning or the use of an ultrasonic cleaner with a short setting of about one minute at a temperature of no more than 40 °C. Cleaning solutions should be added in low concentrations.

Fixing: In principle, any reliable temporary cement can be used, but eugenol-free fixing materials are preferred.

Secondary effects: Such as allergies to contents of the thermoplastic acrylics based reactions may very rarely occur.

Reciprocal actions: unknown

Reactions: In case of known incompatibilities and allergies to contents of the thermoplastic acrylics.

Further questions: Please contact SILADENT technical services (Tel.: +49 (0) 53 21 – 37 79 25 / 26) or our sales representatives.

Polya BioStar PP is classified as a class IIa medical product according to the European directive 93/42/EWG. SILADENT Dr. Böhme & Schöps GmbH packaging carries the CE 0124 mark.

Date of Information: 05/2011

Printing Date: January 15, 2010